
frontier

implante dental de conexión interna hexagonal
dental implant of hexagonal internal connection

www.ilerimplant.com


Made in Spain

frontier

implante dental de conexión interna hexagonal
dental implant of hexagonal internal connection

Sistema de
implantes
Implant system

Índice - Summary

SISTEMA DE IMPLANTES

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
gmi frontier internal connection hexagonal implants are machined in CP grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer that forms on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the rounded shapes, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

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


To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made up of titanium and oxygen, which are the basic components of the oxide layer that forms in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of GMI implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.

SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: Ultrastructural details of cells in biofilm on the implant.

Protocolo de inserción de implantes gmi frontier - Insertion sequence for gmi frontier implants

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gmi frontier internal connection implants come in a screen-printed cardboard box, with the instructions for use, traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi frontier** implants consists of a titanium support which only houses the dental implant, prepared to be inserted directly into the mouth, using a single hex-2.45 mm implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.

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1. Check sticker indicating implant diameter and length and open the tab on the carton.
2. Remove the secondary packaging, instructions for use and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract and save the primary packaging plug making a lateral movement.
6. Insert the HEX-2.45 mm implant carrier wrench until a slight retention is noticed, addressing the notches on the key with the faces of the implant hex broaching.
7. Check that the key is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a maximum torque of 35 N-cm until the treated area is at the crestal or subcrestal level. Remove the key.
9. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
10. Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove key and suture the incision.

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

The following is a description of the drilling sequences for the different internal connection implant models, as well as the recommended drilling conditions for their use:

- Lance-shaped: 1200 - 1500 rpm.
- Pilot drills: 900 - 1200 rpm.
- Guide drill: 800 rpm.
- Twist drills:
 - ∅2.80mm → 500 - 700 rpm.
 - ∅3.00 & ∅3.50mm → 400 - 700 rpm.
 - ∅4.00 & ∅5.10mm → 400 - 600 rpm.
- Countersink drills: 200 - 400 rpm.
- For drilling, external cooling with saline solution is required.


frontier ∅3,30		
	Descr.	Ref.
1	Fresa lanceolada Lance-shaped drill	KYF 0C 1224
2	Fresa piloto Pilot drill	KYF 0C 2221
3	Fresa guía Guide drill	KYF 0C 3227
4	Fresa ∅2,80 mm Twist drill ∅2.80 mm	KYF 0C 0128
5	Fresa ∅3,00 mm Twist drill ∅3.00 mm	KYF 0C 0130
6	Avellanadora ∅3,30 mm Countersink drill ∅3.30 mm	KYF 0C 4133

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



The **gmi frontier** implant range was designed and tested to be inserted as a single restoration, with a 30-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram.

	$\varnothing 3,30\text{-}\varnothing 3,75\text{ mm}$	$\varnothing 4,25\text{ mm}$	$\varnothing 4,75\text{-}\varnothing 5,75\text{ mm}$
	Incisivos (Excepto centrales maxilares) Incisors (Except central maxillary incisors)	Caninos • Canines Premolares • Premolars Incisivos centrales maxilares Central maxillary incisors	Molares Molars

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
The GMI FRONTIER implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. GMI recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The GMI FRONTIER range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore the stability of the connection. GMI recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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
To prevent prior deformation of the clinic screw, GMI recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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Once the repair phase of the tissue supporting the implant has finished, in the case of implants fully or partially covered by soft tissue there must be a mucosal channel or tunnel enabling connection of the implant to the prosthesis or secondary structure. Healing abutments are used to form the mucosal route, once the soft tissue has been opened by incision and the cover screw has been withdrawn, by screwing them to the implant.

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Key features

- Grade V titanium machining as per ISO-5832 standard.
- Tightening with b/f 1.20 mm hex wrench.
- Alternative tightening with a flat tip screwdriver.
- Available in emergency heights from 2.00 to 6.00 mm.
- Apical bit to facilitate the start of threading.

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When the use of a custom abutment is required, millable abutments can be used. By cutting them in any direction, they can be suited to the needs of the prosthetic. In addition and as a complementary function, millable implant abutments can be used like impression coping, threading them to the implant screw through the short impression screw and using them to make the negative of the patient's mouth for subsequent positivation on plaster.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- They allow cutting in any direction to suit the needs of the prosthetic.
- Compatible with b/f 4.00 mm hex wrench for use as a implant carrier.
- Available only in non-rotary version.
- Possibility to use as impression coping with short impression screw.
- Clinic screw part no.: KDT 0C 3601 included.

Plataforma Platform	Ref.
PS · RP	KDB 0F 3601
PA · WP	KDB 0F 3602

Pilares inclinados estéticos gmi frontier - gmi frontier esthetic angled abutments

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
When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used, that once mounted and oriented on the implant will allow the parallelization of the new crown with other structures.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 20 degrees.
- Clinic screw part no.: KDT OC 3601 included.
- Abutment shoulder adapted to the shape of the gum.

Pilares inclinados estéticos Esthetic angled abutments		
Plataforma Platform	Ang. (°)	Ref.
PS · RP	15°	KDM OF 3607
	20°	KDM OF 3611
PA · WP	15°	KDM OF 3609
	20°	KDM OF 3612

Pilares inclinados Angled abutments		
Plataforma Platform	Ang. (°)	Ref.
PS · RP	20°	KDM OF 3613
PA · WP	20°	KDM OF 3614

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Once the mucosal tunnel has been formed by the healing abutment and after removing the cap from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Abutments in b/f 2.45 mm hex non-rotary version.
- Clinic screw part no.: KDT OC 3601 included.
- Available in three emergency heights from 1.50 to 3.50 mm.

Plataforma Platform	Altura (mm) Height (mm)	Ref.
	1,50	KDM OF 3601
PS · RP	2,50	KDM OF 3602
	3,50	KDM OF 3603


Plataforma Platform	Altura (mm) Height (mm)	Ref.
	1,50	KDM OF 3604
PA · WP	2,50	KDM OF 3605
	3,50	KDM OF 3606

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Once the mucosal tunnel has been formed by the healing abutment and after removing the cap from this location, it is appropriate to use the titanium base, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for a metallic or ceramic meso-structure casted or machined with CAD/CAM.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Non-rotary and rotary version.
- Clinic screw part no.: KDT OC 3601 included.
- Available in emergency height of 0.60 mm
- Non-rotary titanium base with indexation.


Plataforma Platform	Tipo conexión Connection type	Ref.
PS • RP	Anti-rotatoria • Non-rotary	KDM OF 3615
	Rotatoria • Rotary	KDM OF 3616
PA • WP	Anti-rotatoria • Non-rotary	KDM OF 3617
	Rotatoria • Rotary	KDM OF 3618

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When it is necessary to perform a provisional restoration while the bone-implant interfaces or the soft tissues are healing, the temporary cylinder is used, defined as technical plastic machined attachment that once attached to the crown of the implant with clinic screw, will absorb most of the masticatory stresses while avoiding overloading the implant while the osseointegration process is completed.

Key features

- Machined in high strength technical plastic.
- Abutment shoulder adapted to the shape of the gum.
- Anti-rotational plane and transverse retention to provide good fixation.
- Clinic screw part no.: KDT OC 3601 included.

Plataforma Platform	Tipo conexión Connection type	Ref.
PS • RP	Anti-rotatoria • Non-rotary	KDH OF 3601
	Rotatoria • Rotary	KDH OF 3603
PA • WP	Anti-rotatoria • Non-rotary	KDH OF 3602
	Rotatoria • Rotary	KDH OF 3604

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
Once the healing abutment has formed the mucosal tunnel and the healing abutment has been removed from this location, the corresponding ball abutment is used, which is defined as a machined attachment being fixed to the top of the implant and serving as the retention element for the removable prosthesis, with the previous addition of the retention caps that will be fixed to the top spherical part of the ball abutment.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.00 mm.
- Pink retention cap and metal container included.
- Protection disk included.
- Tightening with b/f 2.45 mm hex wrench, part no.: KYL 0C 0065.
- Instruments for caps: Universal handle (RH124(CP)), Insertion tool for caps (RH185(IAC)) and Removal tool for caps (RH091(EC)).

Plataforma Platform	Altura (mm) Height (mm)	Ref.
PS - RP	1,00	KDF 0F 3609
	2,00	KDF 0F 3610
	3,00	KDF 0F 3611
	4,00	KDF 0F 3612


Retenciones Retention caps			
Tipo retención Retainer type	Color	Cantidad (ud.) Quantity (Un)	Ref.
Suave - Soft	Rosa - Pink	6	RH40CRNSN
Fuerte - Strong	Verde - Green	6	RH040CRN
Elástica - Elastic	Naranja - Orange	6	RH049PCN

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Following the first stage of surgery and therefore after implant osseointegration is achieved, the multi-esthetic abutment (MEA) is used, that is defined as the machined attachment that, after being attached to the crown of the implant, serves as a support for a multiple restoration, thus decreasing the height between the implant connection and the structure connection, and providing a passive adjustment of the structure.


Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.00 mm.
- One-piece straight abutments in rotary model.
- Tightened with hex keys:
 - RP: B/F 2.00 mm (KYL 0C 0097)
 - WP: B/F 2.50 mm (KYL 0C 0098)
- Impression coping with screw included.
- Castable with clinic screw included.
- Temporary cylinder in grade V titanium with clinic screw included.
- Angled multi-esthetic abutments allow non-parallelism up to 30°.

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


Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

Key features

- Machined replicas in high quality stainless steel.
- Anti-rotation millings and vertical retention neck for a good fixation.

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
Impression coping are defined as the machined attachments that, connecting to the coronal part of the implant within the oral cavity and axially fixed by the impression screw, is used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

Key features

- Machined open tray impression coping in grade V titanium according to ISO-5832.
- Impression screws in grade V titanium according to ISO-5832.
- Tightening of impression screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading impression screw.
- Anti-rotation millings and vertical retention neck for a good fixation.
- Impression coping screw included.
- Long impression screw with longitudinal millings for manual tightening.
- Closed tray impression coping machined in technical plastic.
- Closed tray impression coping whole formed by CT impression coping + hexagonal abutment + screw.

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Tomas impresión cubeta cerrada (CC) Closed tray (CT) impression coping

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
To replace a dental crown with a restoration dental prosthesis, a castable is used, which is defined as a meltable plastic machined attachment that, once positioned over the coronal part of the internal connection implant, will be the initial basis to model the form of the final structure, to be cast in the corresponding metal.

Key features

- Machined castables in white meltable plastic material.
- Physiologically inert and excellent dimensional stability.
- Available in rotary and non-rotary version.
- Clinic screw part no.: KDT OC 3601 included.

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Clinic screws are defined as mechanical devices used for the final anchoring of the abutments, the millable implant abutment or the castable directly to the implant, while laboratory screws are the mechanical devices used for the temporary anchoring of castables to the replica of the implant used by the dental lab where the construction of dental crowns or prostheses takes place.

Key features

- Machined grade V titanium clinic screws according to ISO-5832.
- Machined laboratory screws in stainless steel quality.
- Tightening or screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading.
- Laboratory screw with groove to differ it from clinic screw.

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gmi frontier surgical kit consists of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of gmi frontier internal connection dental implant system. The components are arranged in such a way as to be easy to find, at all times, the necessary instrument at any time for each stage of the intervention.

gmi surgical kits are distributed in cartons, made of high strength technical polymer, and should be steam autoclaved before use at a maximum temperature of 134°C.

KYC OF 3006


	Descr.	Un.	Ref.
Fresa lanceolada	Lance-shaped drill	1	KYF 0C 1224
Fresa piloto ø2,00 mm	Pilot drill ø2.00 mm	1	KYF 0C 2221
Fresa guía ø2,00 / ø2,70 mm	Guide drill ø2.00 / ø2.70 mm	1	KYF 0C 3227
Fresa final ø2,80 mm	Twist drill ø2.80 mm	1	KYF 0C 0128
Fresa final ø3,00 mm	Twist drill ø3.00 mm	1	KYF 0C 0130
Fresa final ø3,50 mm	Twist drill ø3.50 mm	1	KYF 0C 0135
Fresa final ø4,00 mm	Twist drill ø4.00 mm	1	KYF 0C 0140
Fresa final ø5,10 mm	Twist drill ø5.10 mm	1	KYF 0C 0151
Fresa avellanadora frontier ø3,30 mm	Frontier countersink drill ø3.30 mm	1	KYF 0C 4133
Fresa avellanadora frontier ø3,75 mm	Frontier countersink drill ø3.75 mm	1	KYF 0C 4138
Fresa avellanadora frontier ø4,25 mm	Frontier countersink drill ø4.25 mm	1	KYF 0C 4142
Fresa avellanadora frontier ø4,75 mm	Frontier countersink drill ø4.75 mm	1	KYF 0C 4146
Fresa avellanadora frontier ø5,75 mm	Frontier countersink drill ø5.75 mm	1	KYF 0C 4160
Llave carraca IP	TI ratchet wrench	1	KYL 0F 0113
Llave-transportador HEX - 2,45 mm IP corta	Short HEX-2.45 mm implant carrier TI wrench	1	KYL 0F 0096
Llave-transportador HEX - 2,45 mm IP larga	Long HEX-2.45 mm implant carrier TI wrench	1	KYL 0F 0111
Llave-transportador HEX - 2,45 mm CA	HEX-2.45 mm implant carrier HP wrench	1	KYL 0F 0091
Llave carraca corta HEX - 1,20 mm	Short HEX-1.20 mm ratchet wrench	1	KYL 0F 0128
Llave carraca larga HEX - 1,20 mm	Long HEX-1.20 mm ratchet wrench	1	KYL 0F 0129
Alargador llaves CA	HP tip extender	1	KYL 0C 0023
Paralelizador ø2,00 / ø2,80 mm	Paralleling pin ø2.00 / ø2.80 mm	3	KYL 0C 0078
Medidor profundidad ø2,80 mm	Depth gauge ø2.80 mm	1	KYL 0C 0094
Medidor profundidad ø3,50 mm	Depth gauge ø3.50 mm	1	KYL 0C 0095

frontier

implante dental de conexión interna hexagonal
dental implant of hexagonal internal connection

Distribuido por:
Distributed by:

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Pol. Ind. El Segre
C/ Enginyer Mies 705-B
25191 Lleida (SPAIN)
Tel. (+34) 973 184 350
Fax (+34) 973 183 278
www.ilerimplant.com